



Bacterial/ Fungal Lysis Buffer

Product No: K.BUF-9103-50ml

This buffer allows crude nucleic acid extraction from pure isolate of bacterial/ fungal sample. Pure isolate can be from agar culture or cell pellet from liquid culture. This buffer requires addition of Proteinase K enzyme (not provided). The extracted crude DNA is suitable for PCR amplification of 16S rRNA full length gene or ITS gene for Sanger sequencing. It is a solution-based extraction method that utilizes high salt to extract crude nucleic acid for PCR applications.

Storage and Stability

Store the product at room temperature (21 °C – 25 °C). The product expiration date is printed on the label.

Materials Supplied by Users

- ✓ Thermo block, set at 56 °C
- ✓ Centrifuge at 14,000 – 16,000 × g
- ✓ Sterile nuclease-free 1.5 mL micro-centrifuge tubes (2x units per sample)
- ✓ Sterile nuclease-free pipette and pipette tips
- ✓ Isopropanol
- ✓ 70% ethanol
- ✓ Proteinase K solution (20 mg/mL)
- ✓ TE buffer

Protocol

Lysis	<ol style="list-style-type: none"> 1. Sample: <ol style="list-style-type: none"> a) Cut agar: Place 0.5 cm x 0.5 cm of bacterial/ fungal cut agar into a new 1.5 mL micro-centrifuge tube. <div style="text-align: center; border: 1px solid black; width: 100px; margin: 10px auto; padding: 5px;">OR</div> <ol style="list-style-type: none"> b) Transfer 1 mL bacterial/fungal culture into a new 1.5 mL micro-centrifuge tube. Centrifuge at 14,000 × g for 1 minute to collect cell pellet. Discard the supernatant. 2. Add 500 µL of Bacterial/ Fungal Lysis Buffer into 1.5 mL micro-centrifuge tube that contains the recommended sample size of bacteria/fungus. <p><i>Tips:</i></p> <ul style="list-style-type: none"> ✓ 1 mL filter pipette tip is recommended to be used to add Lysis Buffer. ✓ Ensure the entire sample is 100% submerged into Lysis Buffer. ✓ Handle each sample one at a time to avoid cross contamination. 3. Add 3 µL of Proteinase K solution. Vortex to mix and spin down briefly.
-------	--



Lysis	<p>4. Incubate at 56 °C according to incubation time below:</p> <p>a) Bacterial : 30 minutes or overnight (optional)</p> <p>b) Fungal : Overnight</p> <p>5. Centrifuge the lysate at 14,000 – 16,000 × <i>g</i> for 10 minutes.</p>
Washing	<p>6. Transfer ~ 500 μL of supernatant to a new 1.5 mL micro-centrifuge tube, which contains 500 μL of isopropanol. Invert the tube several times to mix gently.</p> <p>7. Centrifuge at 14,000 – 16,000 × <i>g</i> for 10 minutes and discard the supernatant.</p> <p>8. Add 1 mL of 70% ethanol. Centrifuge again at 14,000 – 16,000 × <i>g</i> for 5 minutes and discard the supernatant. <i>Tip: 70% ethanol should be prepared freshly or less than a week.</i></p> <p>9. Air dry the pellet for 3 minutes.</p>
Elution	<p>10. Re-suspend the dried DNA pellet with 50 μL TE Buffer and incubate at 56 °C for < 1 hour. <i>Tip: If necessary, increase the elution volume using not more than 100 μL TE Buffer to dissolve the DNA pellet completely.</i></p> <p>11. Measure the Optical Density (OD) reading using spectrophotometer. Dilute the nucleic acid to 15 – 25 ng/μL. Use 2 μL of this diluted nucleic acid as DNA Template for PCR.</p>

Please contact us at <https://base-asia.com/contact/> for more information.