

Sample Submission Guideline for Standard/ PCR Grade Nucleic Acid Extraction : MBS-6000 series

Sample Type	Standard/ PCR Grade DNA extraction	Standard/ PCR Grade RNA extraction
Tissue/ Cell pellet from Animal, Insect, Bacteria, Fungus, Yeast, Microalgae	<ul style="list-style-type: none"> Amount: 50 - 100 mg No. of replicates: 3 Shipping condition: Dry ice 	<ul style="list-style-type: none"> Amount: 50 - 100 mg No. of replicates: 3 Shipping condition: Dry ice
Plant/ Macroalgae (Seaweed)	<ul style="list-style-type: none"> Amount: 150 - 500 mg wet weight or liquid nitrogen ground powder. Snap freeze fresh plant sample in liquid nitrogen No. of replicates: 3 Shipping condition: Dry ice 	<ul style="list-style-type: none"> Amount: 150 - 500 mg wet weight or liquid nitrogen ground powder. Snap freeze fresh plant sample in liquid nitrogen No. of replicates: 3 Shipping condition: Dry ice
Blood	<p>Option 1: Whole Blood, *Malaysia Customer only</p> <ul style="list-style-type: none"> Volume: up to 3 ml freshly collected blood in EDTA Tube Blood samples to be submitted immediately after withdrawal. Results will not be guaranteed for blood samples that are more than 3 days old No. of replicates: 3 Shipping condition: <ol style="list-style-type: none"> 4°C or Blue ice (reach Apical Scientific within 3 days) Dry ice (reach Apical Scientific > 3 days) <p>Option 2: Buffy Coat</p> <ul style="list-style-type: none"> Collect up to 3 ml blood in EDTA tube, centrifuge the blood with 2,500 x g for 10 minutes using swing bucket centrifuge Remove the plasma with 1 ml syringe, 21G needle Collect the buffy coat in 1.5/2 ml tube No. of replicates: 3 Store buffy coat in -80 °C freezer Shipping condition: Dry ice <p>Option 3: White Blood Cell Pellet</p> <ul style="list-style-type: none"> Volume: up to 3 ml freshly collected blood in EDTA Tube Pretreat the whole blood sample with RBC lysis buffer (1st BASE, K.BUF-9101-100ml), follow Appendix Protocol below for pretreatment. No. of replicates: 3 Store the treated dry WBC pellet in -80 °C freezer. Alternatively, the WBC pellet can be stored in RNALater solution. Shipping condition: Dry ice <p>Appendix Protocol: Pretreatment of whole blood sample with RBC lysis buffer (1st BASE, K.BUF-9101-100ml)</p> <ol style="list-style-type: none"> Transfer 1 mL whole blood into a new 2 mL microcentrifuge tube. Add 1 mL RBC Lysis Buffer (not provided). Invert the tube 10 times. Centrifuge at 700 x g for 5 minutes at room temperature. Carefully remove 1 mL top layer of the supernatant by pipetting. Note: Do not remove the middle and bottom layer which is the white blood cells and red blood cells respectively. Add 1 mL RBC Lysis Buffer (not provided). Resuspend the pellet by pipetting 4 – 5 times. Centrifuge at 700 x g for 5 minutes. Carefully remove 1 mL supernatant by aspirate from top via pipetting. Leave the remaining supernatant and cell pellet in the tube. Repeat Step iv. Remove supernatant completely. [For preservation] Resuspend cell pellet in 1 mL RNALater Solution (not provided). 	<p>White Blood Cell Pellet</p> <ul style="list-style-type: none"> Collect up to 3 ml blood in EDTA tube Must preprocess the blood sample within the same day of sample collection. Pretreat the whole blood sample with RBC lysis buffer (1st BASE, K.BUF-9101-100ml), follow Appendix Protocol below for pretreatment. No. of replicates: 3 Store the treated dry WBC pellet in -80 °C freezer. Alternatively, the WBC pellet can be stored in RNALater solution. Shipping condition: Dry ice <p>Appendix Protocol: Pretreatment of whole blood sample with RBC lysis buffer (1st BASE, K.BUF-9101-100ml)</p> <ol style="list-style-type: none"> Transfer 3 mL whole blood into a 15 mL centrifuge tube. Add 3 mL RBC Lysis Buffer (not provided). Invert the tube 10 times. Centrifuge at 700 x g for 5 minutes at room temperature. Carefully remove 2 mL top layer of the supernatant by pipetting. Note: Do not remove the middle and bottom layer which is the white blood cells and red blood cells respectively. Add 3 mL RBC Lysis Buffer (not provided). Resuspend the pellet by pipetting 4 – 5 times. Centrifuge at 700 x g for 5 minutes. Carefully remove 3 mL supernatant by aspirate from top via pipetting. Leave the remaining supernatant and cell pellet in the tube. Repeat Step iv. Remove supernatant completely. [For preservation] Resuspend cell pellet in 1 mL RNALater Solution (not provided).
Buccal swab sample	<p>Option 1:</p> <ul style="list-style-type: none"> Collect sample using commercial swab sample collection kit Follow collection kit guidelines to preserve the swab sample No. of replicates: 2 Storage condition: Follow the guidelines from the commercial kit Shipping condition: Follow the guidelines from the commercial kit <p>Option 2:</p> <ul style="list-style-type: none"> Collect sample using sterile swab (Optional) Store swab in 1X PBS No. of replicates: 3 Storage condition: -80°C freezer Shipping condition: Dry ice 	<ul style="list-style-type: none"> Collect sample using commercial swab sample collection kit (for RNA application) Follow collection kit guidelines to preserve the swab sample No. of replicates: 2 Storage condition: Follow the guidelines from the commercial kit Shipping condition: Follow the guidelines from the commercial kit
Saliva	<p>Option 1:</p> <ul style="list-style-type: none"> Suggest to follow commercially available saliva collection kit guidelines for sample collection, preservation, storage and shipping. For example, DNA Genotek, OMNIgene, SALIVA DNA and RNA device. No. of replicates: 2 Storage/Shipping condition: Follow kit recommendation <p>Option 2:</p> <ul style="list-style-type: none"> Collect 2 ml saliva in a new 50 ml centrifuge tube Add 10 ml 1X PBS. Shake vigorously to mix for at least 20 seconds. Immediately centrifuge at 2,000 x g for 5 minutes at room temperature to pellet cells. Decant the supernatant without delay. No. of replicates: 2 Store dry pellet in -80 °C freezer Shipping condition: Dry ice 	<p>Option 1:</p> <ul style="list-style-type: none"> Suggest to follow commercially available saliva collection kit guidelines for sample collection, preservation, storage and shipping. For example, DNA Genotek, OMNIgene, SALIVA DNA and RNA device. No. of replicates: 2 Storage/Shipping condition: Follow kit recommendation <p>Option 2:</p> <ul style="list-style-type: none"> Collect 2 ml fresh saliva in a 50 ml centrifuge tube No. of replicates: 2 Add 2 ml RNALater solution. Vortex vigorously to mix thoroughly. Incubate sample at 4°C fridge for overnight. Store in -80°C freezer the next day. Shipping condition: Dry ice
Others	Please enquire.	