

## ExactPro Regular Range (9-180 kDa) Prestained Protein Ladder

**Cat. No.:** BIO-5151-25ul  
BIO-5151-500ul

**Applications:** Used in SDS-PAGE and Western blots for approximation of protein sizes, and to locate a protein of interest for excision from an unstained preparative gel.


**Quality Control:** 5µl of ExactPro Regular Range (9-180 kDa) Prestained Protein Ladder provides 10 blue bands in 4 – 20% gradient SDS-PAGE (Tris-Glycine buffer) and after electro transfer onto nitrocellulose membrane. Two enhanced bands of 25 kDa and 72 kDa observed.

**Storage Buffer:** 20mM Tris-H<sub>3</sub>PO<sub>4</sub> (pH 7.5 @ 25°C)  
15% (v/v) Glycerol  
2% (w/v) SDS  
0.2mM DTT  
3.6mM Urea

**Storage Condition:** 25°C – 2 weeks  
4°C – 3 months  
-20°C – long term storage

### 1.0 DESCRIPTION

ExactPro Regular Range (9-180 kDa) Prestained Protein Ladder contains 10 prestained recombinant proteins that cover a range of molecular weights from 9 kDa – 180 kDa. Proteins are covalently coupled with a blue chromophore, with 2 enhanced reference bands, 25 kDa and 75 kDa, when separated on SDS-PAGE (Tris-Glycine). Prestained proteins may have different mobilities in various SDS-PAGE buffer system.

Band	Color	Tris-Glycine Migration Pattern and Approximate MWs (kDa)	Bis-Tris (MOPS) Migration Pattern and Approximate MWs (kDa)	Bis-Tris (MES) Migration Pattern and Approximate MWs (kDa)	
1	Blue	180	170	170	180
2	Blue	140	130	130	140
3	Blue	100	93	93	100
4	<b>Blue</b>	<b>72</b>	<b>68</b>	<b>70</b>	72
5	Blue	60	53	53	60
6	Blue	45	41	42	45
7	Blue	35	30	30	35
8	<b>Blue</b>	<b>25</b>	<b>22</b>	<b>23</b>	25
9	Blue	15	14	14	15
10	Blue	10	9	9	10

The ladder is ready-to-use. No further dilution, addition of a reducing agent or heating is required.

## 2.0 PROTOCOL

- Thaw ladder at room temperature to dissolve precipitated solids.  
Do not boil.
- Mix solution gently to ensure that it is homogeneous.
- Load ladder into SDS-PAGE gel or Western blot using these volumes (based on gel thickness of 0.75 – 1.0mm).
  - 3  $\mu$ l per loading on 15-well mini-gel
  - 5  $\mu$ l per loading on 10-well mini-gel
  - 1.5 – 2.5  $\mu$ l per loading for general Western transferring
  
- Double the loading volume for thicker or larger gel.