

## REDiant II Taq DNA Polymerase (recombinant)

**Cat. No.** BIO-5116-50U  
BIO-5116-500U  
BIO-5116-2500U

**Concentration:** 0.3U/ $\mu$ l

**Supplied with:** 10X Taq II buffer  
25mM MgCl<sub>2</sub>

**Storage:** -20°C.  
Avoid frequent thawing and freezing.

### Quality Control

**Nuclease Assay:** No detectable contaminating endonuclease or exonuclease activity.

**Functional Assay:** *REDiant II Taq DNA Polymerase was tested for amplification of a 3500bp of human genomic DNA.*

### 1.0 DESCRIPTION

REDiant II Taq DNA Polymerase is a thermostable DNA polymerase that catalyzes a 5'→3' polymerase activity and a 5' flap endonuclease activity. In addition, REDiant II Taq DNA Polymerase exhibits deoxynucleotidyl transferase activity, which frequently results in the addition of extra adenines at the 3'-end of PCR products. Recombinant REDiant II Taq DNA Polymerase is ideal for standard PCR templates up to 6kb. It includes a red, inert tracking dye which migrates at the same rate as 1kb DNA fragment in a 1% agarose gel. Inert dye does not inhibit PCR.

### Source

An *E.coli* strain that carries a cloned gene from *Thermus aquaticus*.

### Unit Definition

One unit of the enzyme catalyzes the incorporation of 10nmol of deoxyribonucleotides into an acid-insoluble material in 30mins at 74°C by using activated salmon sperm DNA as a template/primer.

## Applications

- DNA sequencing
- DNA labeling
- PCR for cloning
- PCR amplification of DNA fragments up to 6kb
- Routine PCR

## Buffer Composition

### Storage buffer

|                    |            |
|--------------------|------------|
| 20mM Tris-HCl      | 0.1mM EDTA |
| Stabiliser         | 1mM DTT    |
| 50% (v/v) Glycerol | 100mM KCl  |

### 10X Taq II buffer

750mM Tris-HCl (pH 8.8 at 25°C)  
200mM KCl  
50mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>  
0.5% (v/v) Nonidet P40

## Remarks

- Half-life of this enzyme is >40mins @ 95°C.
- Taq II DNA Polymerase accepts modified nucleotides (e.g. biotin-, digoxigenin-, fluorescent-labeled nucleotides) as substrates for the DNA synthesis.

## 2.0 PROTOCOL

The following protocol serves as a starting point and general guideline for any Polymerase Chain Reaction (PCR). Reaction conditions (incubation times and temperatures, concentration of Taq II DNA polymerase, primers, MgCl<sub>2</sub> and template DNA) vary and need to be optimised. PCR reactions should be prepared in a DNA-free environment, dedicated pipette and aerosol resistant tips are recommended. Always keep the control and template DNA to be amplified isolated from other components.

**Recommendations of Template DNA in a 50 µl reaction volume**

| DNA Type                  | Amount of DNA |
|---------------------------|---------------|
| Human genomic DNA         | 0.1 to 1 µg   |
| Plasmid DNA               | 0.5 to 5 ng   |
| Phage DNA                 | 0.1 to 10 ng  |
| <i>E.coli</i> genomic DNA | 10 to 100 ng  |

1. Thaw REDiant II Taq DNA polymerase, 10X Taq II buffer and 25mM MgCl<sub>2</sub> at room temperature. Keep the tubes on ice after thawing. Vortex briefly and spin down contents quickly.
2. Prepare the following reaction mix in a sterile, nuclease-free PCR tube on ice.

**Table 1. For a 50µl reaction volume:**

| Components                    | Volume (µl)  | Final Concentration |
|-------------------------------|--------------|---------------------|
| REDiant II Taq DNA Polymerase | 4.5          | 1.5U                |
| Forward Primer (10µM)         | 0.5 – 5      | 0.1 – 1µM           |
| Reverse Primer (10µM)         | 0.5 – 5      | 0.1 – 1µM           |
| 10X Taq II Buffer             | 5            | 1X                  |
| dNTP Mix (10mM of each)       | 1            | 0.2mM               |
| 25mM MgC <sub>2</sub>         | 2-6          | 1 - 3mM             |
| DNA Template                  | 1-5          | See above           |
| Nuclease-Free Water           | Top up to 50 | N.A                 |

3. Cap tubes and spin down contents briefly.
4. Place reactions in thermal cycler and incubate at 95°C to completely denature template DNA. Perform 25 - 35 cycles of PCR amplification by repeating steps 2 to 4 in the PCR reaction. Low amounts of starting template may require 40 cycles.

**Table 2. PCR Amplification as follows:**

| Step                 | Temperature | Time            |
|----------------------|-------------|-----------------|
| Initial Denaturation | 95°C        | 1 – 5 min       |
| Denaturation         | 95°C        | 0.5 – 1 min     |
| Annealing            | 50- 68°C    | 0.5 – 1 min     |
| Extension            | 72 - 75°C   | *0.5 – 1 min/kb |
| Final Extension      | 72 - 75°C   | 5 – 15 min      |
| Soak                 | 4°C         | ∞               |

\* For amplification of >3kb, it is recommended to use 1min/kb as extension rate.

### 3.0 General guidelines for PCR amplifications

#### Initial denaturation

To ensure efficient utilization of the template during the first amplification cycle, it is essential that the template is denatured completely. If GC content of the template is ≤ 50%, an initial 1-5 mins denaturation at 95°C would be sufficient.

### Denaturation

DNA denaturation time of 0.5min/cycle at 95°C is usually sufficient. For GC-rich templates, denaturation could be prolonged to 3-4mins.

### Annealing

Annealing temperature should be 5°C lower than melting temperature ( $T_m - 5^\circ\text{C}$ ) of primers. 0.5min/cycle is usually sufficient. If non-specific PCR products are observed, the temperature should be optimized stepwise in 1-2°C increments.

### Extension

Optimal temperature for extension of Taq II DNA Polymerase is between 70-75°C. Recommended extension step is 30sec/kb at 72°C for PCR products.

### Number of cycles

If less than 10 copies of template are present, about 40 cycles are required. For a higher amount, 25-35 cycles are sufficient.

### Final extension

After the last cycle, incubate PCR mixture at 72°C for an additional 5-15mins to allow any possible incomplete amplification to take place.