

SUPPLEMENTARY PROTOCOL

Ver 1.0

KIT-9022

PrimeWay Genomic II DNA Extraction Kit

A support protocol is an extra guide provided with the main protocol to help users use the product with special or less common sample types, workflows, or conditions. It helps expand the use of the main kit by giving tested methods for other applications, while still ensuring good performance and quality. It doesn't replace the main protocol but offers extra help for different situations.





List of Supplementary Protocol for KIT-9022

1. Formalin-Fixed paraffin Embedded (FFPE) Tissue
2. Formalin-Fixed Tissue
3. Microbial bacteria DNA from saliva



1. Formalin-Fixed paraffin Embedded (FFPE) Tissue

Materials Supplied by User

- ✓ Undenatured absolute ethanol ($\geq 99.5\%$)
- ✓ Nuclease-free water
- ✓ Vortex mixer
- ✓ Centrifuge, at speed of 16,000 – 18,000 x *g*
- ✓ Water bath or dry bath
- ✓ Pipettes & pipette tips
- ✓ 1.5 mL & 2 mL microcentrifuge tubes
- ✓ Xylene

| | |
|-------------|--|
| Preparation | <ol style="list-style-type: none"> I. Set water bath/ dry bath to 37 °C, 60 °C and 90 °C. II. Preheat the Elution Buffer at 60 °C. |
| Sample | <ol style="list-style-type: none"> 1. Transfer up to 25 mg FFPE Tissue to a new 1.5 mL microcentrifuge tube. 2. Add 1.2 mL xylene (not provided) and vortex vigorously. 3. Incubate at room temperature for 10 minutes. 4. Centrifuge at max speed (16,000 – 18,000 x <i>g</i>) for 5 minutes. Remove the supernatant. 5. Add 1.2 mL absolute ethanol (not provided) and gently invert to mix. Centrifuge at max speed (16,000 – 18,000 x <i>g</i>) for 5 minutes. |



| | |
|---------|--|
| Sample | <p>6. Remove the supernatant. Repeat Step 5.</p> <p>7. Remove supernatant completely. Incubate at 37 °C for 10 – 15 minutes, until the excess ethanol has evaporated.</p> |
| Lysis | <p>8. Add 200 µL GL1 Buffer and 20 µL Proteinase K Solution. Vortex to mix.</p> <p>9. Incubate at 60 °C for 1 hour. Invert the tube occasionally.</p> <p>10. Incubate at 90 °C for 1 hour.</p> <p>11. Centrifuge at 14,000 x g for 2 minutes to pellet insoluble debris.</p> <p>12. Transfer the supernatant to a new 1.5 mL microcentrifuge tube.</p> <p>13. Add 200 µL GL2 Buffer. Vortex to mix.</p> |
| Binding | <p>14. Add 200 µL absolute ethanol (not provided). Vortex to mix immediately.</p> <p>15. Place a PrimeWay Genomic II Column into a new Collection Tube.</p> <p>16. Transfer up to 750 µL lysate, including the precipitate if any, to the PrimeWay Genomic II Column. Centrifuge at 14,000 x g for 1 minute.</p> <p>17. Discard the Collection Tube and place the column into a new Collection Tube.</p> |



| | |
|---------|---|
| Washing | <p>18. Add 400 μL Wash Buffer G1 to the column. Centrifuge at 14,000 x <i>g</i> for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> <p>19. Add 600 μL Wash Buffer G2 to the column. Centrifuge at 14,000 x <i>g</i> for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> |
| Drying | <p>20. Centrifuge the column at 14,000 x <i>g</i> for 3 minutes to dry the membrane.</p> |
| Elution | <p>21. Transfer the PrimeWay Genomic II Column to a new 1.5 mL microcentrifuge tube.</p> <p>22. Add 100 μL preheated Elution Buffer to the center of the column membrane. Incubate at room temperature for at least 3 minutes.</p> <p>23. Centrifuge at 14,000 x <i>g</i> for 30 seconds to elute the DNA.</p> |



2. Formalin-Fixed Tissue

Materials Supplied by User

- ✓ Undenatured absolute ethanol ($\geq 99.5\%$)
- ✓ Nuclease-free water
- ✓ Vortex mixer
- ✓ Centrifuge, at speed of 14,000 x g
- ✓ Water bath or dry bath
- ✓ Pipettes & pipette tips
- ✓ 1.5 mL & 2 mL microcentrifuge tubes
- ✓ 1X Phosphate Buffered Saline (PBS)

| | |
|-------------|--|
| Preparation | <ol style="list-style-type: none"> I. Set water bath/ dry bath to 60 °C and 90 °C. II. Preheat the Elution Buffer at 60 °C. |
| Sample | <ol style="list-style-type: none"> 1. Wash the tissue with 1 mL PBS with total of 2 times. Remove the liquid completely. 2. Mince the tissue into tiny pieces with sterile blade and transfer up to 25 mg tissue into a new 1.5 mL microcentrifuge. |
| Lysis | <ol style="list-style-type: none"> 3. Add 200 μL GL1 Buffer and 20 μL Proteinase K Solution. Vortex to mix. 4. Incubate at 60 °C for 1 hour. Invert the tube occasionally. 5. Incubate at 90 °C for 1 hour. |



| | |
|---------|--|
| Lysis | <p>6. Centrifuge at 14,000 x <i>g</i> for 2 minutes to pellet insoluble debris.</p> <p>7. Transfer the supernatant to a new 1.5 mL microcentrifuge tube.</p> <p>8. Add 200 μL GL2 Buffer. Vortex to mix.</p> |
| Binding | <p>9. Add 200 μL absolute ethanol (not provided). Vortex to mix immediately.</p> <p>10. Place a PrimeWay Genomic II Column into a new Collection Tube.</p> <p>11. Transfer up to 750 μL lysate, including the precipitate if any, to the PrimeWay Genomic II Column. Centrifuge at 14,000 x <i>g</i> for 1 minute.</p> <p>12. Discard the Collection Tube and place the column into a new Collection Tube.</p> |
| Washing | <p>13. Add 400 μL Wash Buffer G1 to the column. Centrifuge at 14,000 x <i>g</i> for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> <p>14. Add 600 μL Wash Buffer G2 to the column. Centrifuge at 14,000 x <i>g</i> for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> |
| Drying | <p>15. Centrifuge the column at 14,000 x <i>g</i> for 3 minutes to dry the membrane.</p> |



Elution

16. Transfer the **PrimeWay Genomic II Column** to a new 1.5 mL microcentrifuge tube.
17. Add **100 μ L preheated Elution Buffer** to the center of the column membrane. Incubate at room temperature for at least 3 minutes.
18. Centrifuge at 14,000 x g for 30 seconds to elute the DNA.



3. Microbial bacteria DNA from saliva sample

Materials Supplied by User

- ✓ 1x PBS
- ✓ Bacteria Pre-Lysis Buffer [20 mM Tris-HCl; 2 mM EDTA; 1% Triton X-100; pH 8.0]
- ✓ Lyophilised lysozyme
- ✓ 50 mM EDTA, pH 8.0

| | |
|-------------|--|
| Preparation | <ol style="list-style-type: none"> I. Set water bath/ dry bath to 60 °C. II. Preheat the Elution Buffer at 60 °C. III. Bacteria Pre-Lysis Buffer Mix preparation: Dissolve 20 mg Lysozyme (not provided) into 1 mL Bacteria Pre-Lysis Buffer (not provided). |
| Sample | <ol style="list-style-type: none"> 1. Collect at least 1 mL saliva in a new 50 mL centrifuge tube. 2. Transfer 1 mL saliva to a new 15 mL centrifuge tube. Add 5 mL PBS (not provided). Vortex to mix vigorously. 3. Centrifuge at 2,000 x <i>g</i> for 5 minutes at room temperature to pellet cells. Immediately, discard the supernatant. 4. Weigh the saliva pellet to ensure it is approximately 20–40 mg. |



| | |
|---------|---|
| Sample | <p>Note: <i>If the pellet exceeds 40 mg, double the volume of all subsequent buffers up to the ethanol addition step. Proceed by loading the entire sample onto a single spin column.</i></p> |
| Lysis | <ol style="list-style-type: none"> 5. Resuspend the pellet with 420 µL 50 mM EDTA (not provided). Transfer the sample to a new 1.5 mL microcentrifuge tube. 6. Add 180 µL freshly prepared Bacteria Pre-Lysis Buffer Mix (not provided, refer to preparation method section). Incubate at 37 °C for 30 – 60 minutes. 7. Centrifuge at 16,000 x g for 2 minutes. Discard supernatant. 8. Resuspend the pellet with 200 µL GL1 Buffer. 9. Add 20 µL Proteinase K Solution. Vortex to mix. 10. Incubate at 60 °C for 60 minutes. Invert the tube occasionally. 11. Add 200 µL GL2 Buffer. Vortex to mix. 12. Add 4 µL RNase A Solution. Vortex to mix 13. Incubate at room temperature for 5 minutes. |
| Binding | <ol style="list-style-type: none"> 14. Add 200 µL of absolute ethanol (not provided). Vortex to mix immediately. 15. Place a PrimeWay Genomic II Column into a new Collection Tube. 16. Transfer up to 750 µL lysate, including the precipitate if any, to the PrimeWay Genomic II Column. Centrifuge at 14,000 x g for 1 minute. |



| | |
|---------|---|
| | 17. Discard the Collection Tube and place the column into a new Collection Tube. |
| Washing | <p>18. Add 400 μL Wash Buffer G1 to the column. Centrifuge at 14,000 x g for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> <p>19. Add 600 μL Wash Buffer G2 to the column. Centrifuge at 14,000 x g for 30 seconds. Discard the flow-through and place the column back into the Collection Tube.</p> |
| Drying | 20. Centrifuge the column at 14,000 x g for 3 minutes to dry the membrane. |
| Elution | <p>21. Transfer the PrimeWay Genomic II Column to a new 1.5 mL microcentrifuge tube.</p> <p>22. Add 50 μL preheated Elution Buffer to the centre of the column membrane. Incubate at room temperature for at least 3 minutes.</p> <p>23. Centrifuge at 14,000 x g for 30 seconds to elute the DNA.</p> |

Product Ordering Information

| Protocol | Part Number | Product Description | Remarks |
|--------------|----------------------------------|--|---------|
| FFPE tissues | 1st BASE BUF-2041- 1x500ml | 1X Phosphate Buffered Saline (PBS), Biotechnology Grade, 500 mL | |