

Product Specification Sheet



Product Code : BUF-1170

Identifier : Water, DEPC treated

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : H₂O

Molecular Weight : 18.02 g/mol

Product Description : Water incubated with Diethylpyrocarbonate (DEPC) and autoclaved to remove DEPC.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com