

Product Specification Sheet



Product Code : BUF-1180

Identifier : Water

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : H₂O

Molecular Weight : 18.02 g/mol

Product Description : Sterilized by 0.1µm membrane filtration, followed by autoclaving. Suitable for molecular biology applications.

Storage Condition : Store at room temperature.

Assay:

Appearance
DNase/ RNase activity

Specification:

Clear Liquid
None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com