

# Product Specification Sheet



**Product Code** : BUF-1415  
**Identifier** : 1.0M Tris Buffer, pH 6.8  
**Grade** : Ultra Pure  
**Appearance** : Clear, colorless liquid  
**Molecular Formula** : Not Applicable  
**Molecular Weight** : Not Applicable  
**Product Description** : 1.0M Tris buffer solution adjusted to pH 6.8.  
**Storage Condition** : Store at room temperature.

<b><u>Assay:</u></b>	<b><u>Specification:</u></b>
Appearance pH (neat, 25°C)	Clear Liquid 6.6 – 7.0

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com