

Product Specification Sheet



Product Code : BUF-1419
Identifier : 1.5M Tris Buffer, pH 8.8
Grade : Biotechnology
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : 1.5M Tris buffer solution adjusted to pH 8.8.
Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	8.6 – 9.0
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com