

# Product Specification Sheet



**Product Code** : BUF-2040

**Identifier** : 1X Phosphate Buffered Saline (PBS)

**Grade** : Ultra Pure

**Appearance** : Clear, colorless liquid

**Molecular Formula** : Not Applicable

**Molecular Weight** : Not Applicable

**Product Description** : Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 12mM Phosphate Buffer at 1X concentration.

**Storage Condition** : Store at room temperature.

**Assay:**

Appearance  
pH (neat, 25°C)

**Specification:**

Clear Liquid  
7.2 – 7.6

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com