

Product Specification Sheet



Product Code : BUF-2041

Identifier : 10X Phosphate Buffered Saline (PBS)

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 12mM Phosphate Buffer at 1X concentration.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	6.5 – 6.9
pH (1X, 25°C)	7.2 – 7.6
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com