

# Product Specification Sheet



**Product Code** : BUF-3001  
**Identifier** : 10X Tris-Acetate-EDTA (TAE) Buffer, pH 7.8  
**Grade** : Ultra Pure  
**Appearance** : Clear, colorless liquid  
**Molecular Formula** : Not Applicable  
**Molecular Weight** : Not Applicable  
**Product Description** : Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.  
**Storage Condition** : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	7.60 – 8.00
Conductivity (1X Conc.)	1100 - 1750 $\mu$ S

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com