

# Product Specification Sheet



**Product Code** : BUF-3020  
**Identifier** : 10X Tris-EDTA (TE) Buffer, pH 8.0  
**Grade** : Ultra Pure  
**Appearance** : Clear, colorless liquid  
**Molecular Formula** : Not Applicable  
**Molecular Weight** : Not Applicable  
**Product Description** : Contains 10mM Tris-HCl and 1mM EDTA at 1X concentration.  
**Storage Condition** : Store at room temperature.

**Assay:**

Appearance  
pH (neat, 25°C)

**Specification:**

Clear Liquid  
7.8 – 8.2

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com