

# Product Specification Sheet



**Product Code** : BUF-3060

**Identifier** : 1X SSPE Buffer

**Grade** : Biotechnology

**Appearance** : Clear, colorless liquid

**Molecular Formula** : Not Applicable

**Molecular Weight** : Not Applicable

**Product Description** : Contains 0.15M Sodium Chloride, 0.01M Sodium Phosphate and 1mM EDTA at 1X concentration. Suitable as washing buffer for microarrays.

**Storage Condition** : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	7.2 – 7.6
DNase/ RNase activity	None Detected

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com