


UNLOCK The DNA Codes of Microbes

Transforming microbial identification from a traditional species-based to gene-based view.

FREE 10 rxns
Sequencing+ PLUS Reactions
for every purchase of a
Bacterial/Fungal DNA Barcoding PCR kit
Limited Time Offer!



1st BASE's DNA Barcoding Kits were designed to amplify conserved region of bacterial 16s rRNA gene and fungal ITS gene respectively for Sanger sequencing. It comes with the option to enable crude nucleic acid extraction using high salt solution-based extraction method within **30 minutes** or overnight incubation.

| Product No. | Product Description | Size |
|-------------|---------------------------------|----------|
| KIT-1100-50 | Bacterial DNA Barcoding Kit | 50 Preps |
| KIT-1110-50 | Bacterial DNA Barcoding PCR Kit | 50 Preps |
| KIT-1200-50 | Fungal DNA Barcoding Kit | 50 Preps |
| KIT-1210-50 | Fungal DNA Barcoding PCR Kit | 50 Preps |



DNA Barcoding Workflow

