

Product Information

Binding Enhancer for Plant Lysis Buffer, 1mL

C/No. K.RGT-9101-1ml
Packaging 1 mL
Storage - 20 °C

Apical Scientific Sdn Bhd
7-1 to 7-4, Jalan SP 2/7
Taman Serdang Perdana,
43300 Seri Kembangan
Selangor, Malaysia



Page 1 of 3

Description

Sterilized with 0.22 μ M membrane filtration. This buffer is added into the Plant Lysis Buffer. It helps in the binding and precipitation of nucleic acids during extraction.

Applications

- Nucleic acid extraction

Page 2 of 3

Quality Control

DNA Purity Assay

Genomic DNA was extracted from plant. Absorbance with the ratio 260/280 between 1.7 – 2.0 was obtained. Single intact genomic DNA band was observed.

PCR Amplification Assay

The *rbcl* region was amplified from the purified plant genomic DNA.

Page 3 of 3