

Product Information

Lithium Chloride Solution, 8 M

C/No. K.RGT-9102-10ml
Packaging 10 mL
Storage Room temperature
(21 °C – 25 °C)

Apical Scientific Sdn Bhd
7-1 to 7-4, Jalan SP 2/7
Taman Serdang Perdana,
43300 Seri Kembangan
Selangor, Malaysia



Page 1 of 3

Description

Lithium Chloride Solution is a ready-to-use and nuclease-free solution for RNA isolation. This solution allows precipitation of RNA but not DNA, protein, or free nucleotides triphosphates.

Applications

- RNA precipitation

Concentration

8 mol/L

Page 2 of 3

Quality Control

DNase Assay

No detectable contaminating endonuclease or exonuclease activity.

RNase Activity

RNase contaminants not detected.

Page 3 of 3