

# Product Specification Sheet



**Product Code:** BUF-1052-pH8.0  
**Identifier:** 0.5M EDTA Solution, pH 8.0  
**Grade:** Ultra Pure Grade  
**Appearance:** Clear, colorless liquid  
**Molecular Formula:** Not Applicable  
**Molecular Weight:** Not Applicable  
**Description:** A pH-adjusted EDTA solution with NaOH.  
**Storage Condition:** Store between 20°C - 25°C  
**Additional Info:** Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	7.80 – 8.20

## Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305  
Tel: +65-6775 7318 Email: [custcare@axilscientific.com](mailto:custcare@axilscientific.com)  
Website: [www.axilscientific.com](http://www.axilscientific.com)