

Product Specification Sheet



Product Code: BUF-1151
Identifier: 3.0M Sodium Acetate Solution, pH 5.2
Grade: Biotechnology Grade
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Not Applicable
Storage Condition: Store between 20°C - 25°C
Additional Info: Not Applicable

| Assay: | Specification: |
|-----------------------|----------------|
| Appearance | Clear Liquid |
| pH (neat, 25°C) | 5.00 – 5.40 |
| DNase/ RNase activity | None Detected |

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com