

Product Specification Sheet



Product Code : BIO-1200
Identifier : Potassium Dihydrogen Phosphate
Grade : ACS
Appearance : White powder
Molecular Formula : KH_2PO_4
Molecular Weight : 136.09 g/mol
Product Description : -
Storage Condition : Store at room temperature. Protect from moisture.

<u>Assay:</u>	<u>Specification:</u>
Purity (KH_2PO_4)	$\geq 99.5 \%$
Solubility in water	PASS
Loss on drying (105°C)	$\leq 0.2 \%$
pH (50g/L, 25°C)	4.2 – 4.5
Chloride (Cl)	$\leq 5 \text{ ppm}$
Nitrate	PASS
Sulfate (SO_4)	$\leq 0.003 \%$
Heavy metals (as Pb)	$\leq 0.001 \%$
Sodium (Na)	$\leq 0.02 \%$
Arsenic (As)	$\leq 1 \text{ ppm}$
Iron (Fe)	$\leq 0.001 \%$
Ammonium (NH_4)	$\leq 0.001 \%$

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com