

Product Specification Sheet



Product Code : BIO-1300
Identifier : Potassium Chloride
Grade : ACS
Appearance : White crystals
Molecular Formula : KCl
Molecular Weight : 74.55 g/mol
Product Description : -
Storage Condition : Store at room temperature. Protect from moisture.

<u>Assay:</u>	<u>Specification:</u>
Purity (KCl)	≥ 99.5 %
Solubility in water	P/F
Loss on Drying (110°C)	≤ 0.1 %
pH (50g/L, 25°C)	5.0 – 8.0
Bromide	PASS
Iodide	PASS
Phosphate (PO ₄)	≤ 5 ppm
Sulfate (SO ₄)	≤ 0.002 %
Nitrogen compounds (as N)	≤ 0.001 %
Sodium (Na)	≤ 0.02 %
Copper (Cu)	≤ 2 ppm
Magnesium (Mg)	≤ 0.002 %
Calcium (Ca)	≤ 0.002 %
Barium (Ba)	≤ 0.001 %
Lead (Pb)	≤ 2 ppm
Iron (Fe)	≤ 2 ppm

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com