

Product Specification Sheet



Product Code : CUS-1416
Identifier : 1.0M Tris Buffer, pH 7.4
Grade : Biotechnology
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : 1.0M Tris buffer solution adjusted to pH 7.4.
Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	7.2 – 7.6
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com