

Product Specification Sheet



Product Code: BUF-1419, pH6.8
Identifier: 1.5M Tris Buffer, pH 6.8
Grade: Ultra Pure Grade
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: 1.5M Tris buffer solution adjusted to pH 6.8.
Storage Condition: Store between 20°C - 25°C
Additional Info: Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	6.60 - 7.00

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com