

# Product Specification Sheet



**Product Code** : BUF-1419  
**Identifier** : 1.5M Tris Buffer, pH 8.8  
**Grade** : Ultra Pure  
**Appearance** : Clear, colorless liquid  
**Molecular Formula** : Not Applicable  
**Molecular Weight** : Not Applicable  
**Product Description** : 1.5M Tris buffer solution adjusted to pH 8.8.  
**Storage Condition** : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance pH (neat, 25°C)	Clear Liquid 8.6 – 9.0

**Axil Scientific Pte Ltd**

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610  
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com  
Website: www.axilscientific.com