

Product Specification Sheet



Product Code: BUF-2000-10X
Identifier: 10X Tris Tricine Buffer
Grade: Ultra Pure Grade
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 0.1M Tris and 0.1M Tricine at 1X concentration.
Storage Condition: Store between 20°C - 25°C
Additional Info: Suitable for gel electrophoresis with Tricine gels.

Assay:	Specification:
Appearance	Clear Liquid
Colour	Colourless

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com