

Product Specification Sheet



Product Code : BUF-2020

Identifier : 1X Tris Glycine (TG) Buffer

Grade : Ultra Pure

Appearance : Clear, pale yellow liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 25mM Tris and 192mM Glycine at 1X concentration. Suitable for Western blotting and gel electrophoresis.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
Color	Pale Yellow
pH (1X, 25°C)	8.3 - 8.7

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com