

Product Specification Sheet



Product Code: BUF-2020-1X
Identifier: 1X Tris Glycine (TG) Buffer
Grade: Ultra Pure Grade
Appearance: Clear, pale yellow liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 25mM Tris and 192mM Glycine at 1X concentration.
Storage Condition: Store between 20°C - 25°C
Additional Info: Suitable for Western blotting and gel electrophoresis.

Assay:	Specification:
Appearance	Clear Liquid
Colour	Pale Yellow
pH (1X, 25°C)	8.30 - 8.70

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com