

# Product Specification Sheet



**Product Code:** BUF-2040-10X  
**Identifier:** 10X Phosphate Buffered Saline (PBS)  
**Grade:** Ultra Pure Grade  
**Appearance:** Clear, colorless liquid  
**Molecular Formula:** Not Applicable  
**Molecular Weight:** Not Applicable  
**Description:** Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 12mM Phosphate Buffer at 1X concentration.  
**Storage Condition:** Store between 20°C - 25°C  
**Additional Info:** Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	6.50 - 6.90
pH (1X, 25°C)	7.20 – 7.60

## Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305  
Tel: +65-6775 7318 Email: [custcare@axilscientific.com](mailto:custcare@axilscientific.com)  
Website: [www.axilscientific.com](http://www.axilscientific.com)