

Product Specification Sheet



Product Code: BUF-2041-10X
Identifier: 10X Phosphate Buffered Saline (PBS)
Grade: Biotechnology Grade
Appearance: Clear, colorless liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 12mM Phosphate Buffer at 1X concentration.
Storage Condition: Store between 20°C - 25°C
Additional Info: Sterilized by autoclaving.

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	6.50 – 6.90
pH (1X, 25°C)	7.20 – 7.60
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com