

Product Specification Sheet



Product Code : BUF-3000

Identifier : 10X Tris-Acetate-EDTA (TAE) Buffer, pH 8.0

Grade : Ultra Pure

Appearance : Clear, colorless liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.

Storage Condition : Store at room temperature.

Assay:

Appearance
pH (neat, 25°C)
Conductivity (1X Conc.)

Specification:

Clear Liquid
7.80 – 8.20
1320 - 1980 μ S/cm

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com