

# Product Specification Sheet



**Product Code** : BUF-3000  
**Identifier** : 50X Tris-Acetate-EDTA (TAE) Buffer, pH 8.0  
**Grade** : Ultra Pure  
**Appearance** : Clear, colorless liquid  
**Molecular Formula** : Not Applicable  
**Molecular Weight** : Not Applicable  
**Product Description** : Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.  
**Storage Condition** : Store at room temperature.

<b><u>Assay:</u></b>	<b><u>Specification:</u></b>
Appearance	Clear Liquid
pH (neat, 25°C)	7.80 – 8.20
Conductivity (1X Conc.)	1320 - 1980 $\mu$ S/cm

**Axil Scientific Pte Ltd**

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305  
Tel: +65-6775 7318 Email: [custcare@axilscientific.com](mailto:custcare@axilscientific.com)  
Website: [www.axilscientific.com](http://www.axilscientific.com)