

Product Specification Sheet



Product Code : BUF-3000
Identifier : 50X Tris-Acetate-EDTA (TAE) Buffer, pH 8.0
Grade : Ultra Pure
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.
Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	7.80 – 8.20
Conductivity (1X Conc.)	1100 - 1750 μ S/cm

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com