

Product Specification Sheet



Product Code : BUF-3003
Identifier : 10X Tris-Acetate-EDTA (TAE) Buffer, pH 8.0
Grade : Biotechnology
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.
Storage Condition : Store at room temperature.

Assay:

Appearance
pH (neat, 25°C)
Conductivity (1X Conc.)
DNase/ RNase activity

Specification:

Clear Liquid
7.80 – 8.20
1100 - 1750µS
None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com