

Product Specification Sheet



Product Code: BUF-3003-10X
Identifier: 10X Tris-Acetate-EDTA (TAE) Buffer, pH 8.0
Grade: Biotechnology Grade
Appearance: Clear, colorless liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 40mM Tris-Acetate and 1mM EDTA at 1X concentration.
Storage Condition: Store between 20°C - 25°C.
Additional Info: Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	7.80 – 8.20
Conductivity (1X Conc.)	1320 - 1980 μ S/cm
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com