

Product Specification Sheet



Product Code : BUF-3021
Identifier : 1X Tris-EDTA (TE) Buffer with reduced EDTA, pH 8.0
Grade : Ultra Pure
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : Contains 10mM Tris-HCl and 0.1mM EDTA at 1X concentration.
Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance pH (neat, 25°C)	Clear Liquid 7.8 – 8.2

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com