

Product Specification Sheet



Product Code : BUF-3024
Identifier : 10X Tris-EDTA (TE) Buffer, pH 7.0
Grade : Biotechnology
Appearance : Clear, colorless liquid
Molecular Formula : Not Applicable
Molecular Weight : Not Applicable
Product Description : Contains 10mM Tris-HCl and 1mM EDTA at 1X concentration.
Storage Condition : Store at room temperature.

Assay:

Appearance
pH (neat, 25°C)
DNase/ RNase activity

Specification:

Clear Liquid
6.8 – 7.2
None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com