

Product Specification Sheet



Product Code : BUF-3030

Identifier : 10X Tris Buffered Saline (TBS) Buffer

Grade : Ultra Pure

Appearance : Clear, colorless liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 0.14M Sodium Chloride, 3mM Potassium Chloride and 25mM Tris Base at 1X concentration.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance pH (neat, 25°C)	Clear Liquid 7.20 – 7.60

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com