

Product Specification Sheet



Product Code: BUF-3030-20X
Identifier: 20X Tris Buffered Saline (TBS) Buffer
Grade: Ultra Pure Grade
Appearance: Clear, colorless liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 0.14M Sodium Chloride, 3mM Potassium Chloride and 25mM Tris Base at 1X concentration.
Storage Condition: Store between 20°C - 25°C.
Additional Info: Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	7.20 – 7.60

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com