

Product Specification Sheet



Product Code : BUF-3050

Identifier : 1X SSC Buffer

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 0.15M Sodium Chloride and 15mM Sodium Citrate at 1X concentration. Commonly used in Southern and Northern blotting hybridization and washing procedures.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	6.8 – 7.2
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

41 Science Park Road, #04-08 The Gemini, Singapore Science Park II, Singapore 117610
Tel: +65-6775 7318 Fax: +65-6775 7211 Email: info@axilscientific.com
Website: www.axilscientific.com