

Product Specification Sheet



Product Code:	BUF-3050-1X
Identifier:	1X SSC Buffer
Grade:	Biotechnology Grade
Appearance:	Clear Liquid
Molecular Formula:	Not Applicable
Molecular Weight:	Not Applicable
Description:	Contains 0.15M Sodium Chloride and 15mM Sodium Citrate at 1X concentration.
Storage Condition:	Store between 20°C - 25°C.
Additional Info:	Commonly used in Southern and Northern blotting hybridization and washing procedures.

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	6.80 – 7.20
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com