

Product Specification Sheet



Product Code: BUF-3050-20X
Identifier: 20X SSC Buffer
Grade: Biotechnology Grade
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Contains 0.15M Sodium Chloride and 15mM Sodium Citrate at 1X concentration.
Storage Condition: Store between 20°C - 25°C.
Additional Info: Commonly used in Southern and Northern blotting hybridization and washing procedures.

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	6.80 – 7.20
DNase/ RNase activity	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com