

# Product Specification Sheet



**Product Code:** BUF-4100  
**Identifier:** Deionized Water  
**Grade:** Biotechnology Grade  
**Appearance:** Clear Liquid  
**Molecular Formula:** H<sub>2</sub>O  
**Molecular Weight:** 18.02g/mol  
**Description:** 0.1µm membrane filtered. Suitable for molecular biology applications.  
**Storage Condition:** Store between 20°C - 25°C.  
**Additional Info:** Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
DNase/ RNase activity	None Detected

## Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305  
Tel: +65-6775 7318 Email: [custcare@axilscientific.com](mailto:custcare@axilscientific.com)  
Website: [www.axilscientific.com](http://www.axilscientific.com)