

# Dual-Labeled Fluorogenic Probes for Quantitative Real-Time PCR

- HPLC purified
- Min. Guaranteed 5nmoles delivered quantity
- All probes are quality checked by mass spectrometry
- 15 - 35 Bases\*
- Delivered in 7 working days

## Available Labels

5' Reporter	3' Quencher
6-FAM	TAMRA
HEX	BHQ-1
TET	BHQ-2
TAMRA	BHQ-3
CY5	DABCYL
CY3	QXL 490
ROX	QXL 570
JOE	QXL 670
TEXAS RED	

\* 35 - 45 mers length also available, conditions apply

LNA probes also available, please enquire

Place your order at  
[oligo-sg@base-asia.com](mailto:oligo-sg@base-asia.com)

For enquiries, please contact

Joanna Leong  
HP: 8123 6383  
[joanna-leong@base-asia.com](mailto:joanna-leong@base-asia.com)

Rachel Wong  
HP: 9456 3027  
[rachel-wong@base-asia.com](mailto:rachel-wong@base-asia.com)

Cortina Muljadi  
HP: 9879 8844  
[cortina@base-asia.com](mailto:cortina@base-asia.com)



1st BASE Pte Ltd  
41 Science Park Road #01-23 The Gemini  
Singapore Science Park II, Singapore 117610  
Tel: +65 6775 7318  
Fax: +65 6775 7211  
Website: [www.base-asia.com](http://www.base-asia.com)