

SPECIAL OFFER
on IDT Oligo
Modifications &
Purifications!

GET 15% off*
Oligo Purification,
Modification OR BOTH

Promo Period: 2 Jul - 28 Sep 2012

Every purified oligo up to 60 bases in length includes an additional QC by Capillary Electrophoresis or HPLC (on top of the standard mass spectrometry) to validate final purity.

Types of IDT Oligo Purifications Available

Options	Suited For	Purity
PAGE	5 - 100 bases	> 90 %
HPLC	5 - 60 bases	> 85 %
RNase-free HPLC	Applications sensitive to Ribonucleases	> 85 %

Types of IDT Oligo Modifications Available

Fluorophores

Alexa Fluor® Dyes, Cy™ Dyes, 6-FAM™, HEX™, TET™, TAMRA™, Fluorescein dT, Texas Red® -X (NHS Ester), Rhodamine Dyes, JOE (NHS Ester), LI-COR IRDyes®, WellRED Dyes, Lightcycler® 640 (NHS Ester) and more

Freedom Dyes

6-FAM™ (Fluorescein), MAX (NHS Ester), TYE™ Dyes, TEX 615

Dark Quenchers

Iowa Black®, Black Hole Quencher®, Dabcyl

Attachment Chemistry / Linkers & Spacers

Acrydite, Amino Modifiers, Biotinylation, Thiol, C3 Spacers, Digoxigenin (NHS Ester), Spacer 9, PC Spacer, Cholesteryl-TEG and more

Modified Bases

2-Aminopurine, DeoxyInosine, DeoxyUridine, Inverted dT, 5' Nitrodole, Trimer-20, Fluorine Bases, 2'-O-Methyl RNA Bases and more

Phosphorylation

Phosphorothioate Bonds

*Offer applies to the listed options of purification and modification, enquire for details

100%

IDT offers 100% QC and Purity Guarantees.

All mass spectrometry traces are available for FREE on the IDT website!

First BASE Laboratories Sdn Bhd is the exclusive distributor for Integrated DNA Technologies (IDT) in Malaysia.

Order online at <http://eu.idtdna.com/order/order.aspx> OR contact us directly.

First BASE Laboratories Sdn Bhd

No 7-1 to 7-3, Jalan SP 2/7, Tmn Serdang Perdana Seksyen 2, Seri Kembangan 43300
Selangor Darul Ehsan, Malaysia

Tel: +603 8943 3252

Fax: +603 8943 3243

Email: oligo-my@base-asia.com