

Perform Your Nucleic Acid Extraction With Confidence!




1st BASE

YOUR RESEARCH, OUR BUSINESS

CONVENIENCE


Simple and straight-forward procedure that saves time and eliminates manual processing steps.



Easy-to-follow 1st BASE PrimeWay spin column purification method.

HIGH YIELD


Ensures dependable data with excellent yield.



DNA yield of PCR product using 1st BASE PrimeWay Gel Extraction/Purification Kit.

PURITY CONSISTENCY



Purity ratio of DNA remains consistent within the acceptable range in multiple samplings.



DNA purity assessment from different sample matrices using 1st BASE PrimeWay Plasmid DNA Extraction Kit.

VERSATILITY

Every PrimeWay extraction kit is able to accommodate to an extensive range of sample types and applications.






Fig.1 - KIT-9040-50

Maximum Efficiency & Performance Throughout Your Extraction Processes!



Find Out More Here!

Extraction Kits	Viral DNA/RNA	Genomic DNA	Plant DNA	Plasmid DNA	Gel Extraction/PCR Purification
Product Code	KIT-9010-50	KIT-9020-50	KIT-9030-50	KIT-9040-50	KIT-9050-50
Size	50 preps	50 preps	50 preps	50 preps	50 preps
Sample Type	Serum Plasma Body Fluids Cell Culture Supernatant Viral Transport Medium Swabs	Tissue/ Mouse Tails Buccal Swab Bacteria Cultured Cells Dried Blood Spots Yeast Insects Blood/ Buffy Coat	Various Types of Plants	Bacteria Culture	Gel/ PCR mixture
Yield/ Recovery	<1 µg	20 - 35 µg	1 - 30 µg	Up to 40 µg	Up to 95%
Purity	-	1.7 - 2.0	1.7 - 2.0	1.7 - 2.0	1.7 - 2.0
Duration	30 mins	20 mins	30 mins	25 mins	Gel Extraction: 20 mins PCR Purification: 10 mins



**Find the Right Extraction Kits for your Research Needs,
Contact Us Today!**



SINGAPORE

+65-67757318

✉ custcare@axilscientific.com

🌐 www.axilscientific.com



MALAYSIA

+603-89433252

✉ custcare@apicalscientific.com

🌐 www.apicalscientific.com