

Product Specification Sheet



Product Code : BUF-2046

Identifier : 1X Dulbecco's Phosphate Buffered Saline (D-PBS), without Ca and Mg

Grade : Biotechnology

Appearance : Clear Liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 9.6mM Phosphate Buffer at 1X concentration. Sterilized by 0.2µm membrane filtration.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
pH (neat, 25°C)	7.2 – 7.6
DNase/ RNase activity	None Detected
Endotoxin Level (EU/ml)	< 0.25

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com