

Product Specification Sheet



Product Code : BUF-2048

Identifier : 1X Phosphate Buffered Saline, pH 7.2, Sterile Filtered

Grade : Biotechnology

Appearance : Clear Liquid

Molecular Formula : Not Applicable

Molecular Weight : Not Applicable

Product Description : Contains 155mM Sodium Chloride, 2.71mM Sodium Phosphate, Dibasic and 1.54mM Potassium Phosphate, Monobasic. Sterilized by 0.2µm membrane filtration and autoclaved.

Storage Condition : Store at room temperature.

Assay:

Appearance
pH (neat, 25°C)
DNase/ RNase activity
Endotoxin Level (EU/ml)

Specification:

Clear Liquid
7.0 – 7.4
None Detected
< 0.25

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com