

Product Specification Sheet



Product Code: BUF-5210-1X
Identifier: 1X RIPA Buffer
Grade: Biotechnology Grade
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Consists of 50 mM Tris pH 7.6, 150 mM NaCl, 1% NP-40, 0.5% sodium deoxycholate, 0.1% SDS at 1X concentration.
Storage Condition: Store dark in 4°C.
Additional Info: Not Applicable

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	7.40 - 7.80
Conductivity (neat, 25°C)	As recorded
DNase/ RNase activity	None Detected

Designed by:	Noreeman (Operations Manager)		Customer Approval (Name, Signature, Date):	
Checked by:	Alice Wong Sin Yee (Quality Manager)			
Date:	10-Jun-2024			

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com