

1st BASE Powered by Nacalai Tesque Cell Culture Media



Animal Origin-Free

Endotoxin Tested

pH 7.1-7.5

Sterility Tested for Bacteria, Fungus and Mycoplasma

Product Name	Product No.	Packaging
DMEM (4.5g/L Glucose) with L-Gln, with Sodium Pyruvate, liquid	CUL-5000	500ml, 1L
DMEM (4.5g/L Glucose) with L-Gln, without Sodium Pyruvate, liquid	CUL-5010	500ml, 1L
DMEM/Ham's F-12 with L-Gln, Sodium Pyruvate and HEPES, without Phenol Red, liquid	CUL-5100	500ml, 1L

*Please indicate your size preference of either 5ml or 10ml when placing your order. Serological pipettes are individually wrapped in paper-plastic, sterile, and come in boxes of 200.



1st BASE is proud to announce its collaboration with Nacalai Tesque, a Japan-based company with over 60 years of expertise in cell culture media.

This partnership combines Nacalai Tesque's rich history and field expertise with the efficient and cost-effective capabilities of 1st BASE Contract Manufacturing (CM). Together, we aim to deliver high-quality cell culture media to Southeast Asia, adhering to 1st BASE CM's stringent manufacturing standards.

