

Product Specification Sheet



Product Code	:	K.BUF-9107
Identifier	:	Stool DNA Collection Kit
Appearance	:	Clear, colorless liquid
Molecular Formula	:	Not Applicable
Molecular Weight	:	Not Applicable
Product Description	:	Stool DNA Collection Kit is designed for collecting, preserving, and storing stool samples for DNA extraction. It contains a non-lysis solution to ensure the stability of sample.
Shipping Condition	:	Room temperature
Storage Condition	:	Room temperature

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear, colorless liquid
pH	11.20– 11.50
Conductivity	7.7 – 8.0 mS/cm

For Research Use Only. Not for use in Diagnostic Procedures.

Apical Scientific Sdn Bhd (200301002524)

Lot 7-1 to 7-4, Jalan SP 2/7, Taman Serdang Perdana, Seksyen 2, Seri Kembangan 43300, Selangor, Malaysia
Tel: +603-8943 3252 Fax: +603-8943 3243 Email: custcare@apicalscientific.com Website: www.base-asia.com